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The screenshot shows the 'Sieving' configuration window in the Retsch software. The 'Sample weight' is set to 3000 g. Under 'Sample tolerance and check', the 'Weigh-in tolerance' is set to 'On' with a range from -5.00 % to 5.00 %, and 'Weighing assistant' is set to 'Off'.

AIR JET SIEVING MACHINES AS 200 JET PRO / AS 200 JET PHARMA

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
## TREND ANALYSIS: TEST EQUIPMENT AND PROCESS MONITORING - PROACTIVELY AVOID INCORRECT SIEVING

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
Trend Analysis with downward outliers of 5 consecutive sieve analyses (blinding)

AIR JET SIEVING MACHINES AS 200 JET PRO / AS 200 JET PHARMA

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
Fast sieving   Manual sieving   **Methods**

Search  + Add method

Name	Sample material	Actions
Coffee Arabica	Ground Coffee	 
Coffee Barista	Ground Coffee	 
LAST USED	Ground Coffee	 

**1/11 OPTIONAL SETUP**  
**CHOOSE A METHOD FOR EDITING**

Sieving

Fast sieving Manual sieving **Methods**

Search Cancel    

Basic **Sample** Methodology Data Device

Define sample information

**Characteristics**

Sample material Ground Coffee

Sample preparation No sample preparation ▾

Density 0.37 g/cm<sup>3</sup> 

Source Coffee Mill 04

Sampling

Sample weight 0.00 g ▾

**Sample tolerances and check**

Off  Weigh-in tolerance -5.00 % ▾ To 5.00 % ▾

**On**  Weighing assistant

Comment

**2/11 OPTIONAL SETUP  
ACTIVATION OF THE WEIGHING ASSISTANT**

Sieving

Fast sieving   Manual sieving   **Methods**

Search  + Add method

Name	Sample material	Actions
Coffee Arabica	Ground Coffee	 
Coffee Barista	Ground Coffee	 
LAST USED	Ground Coffee	 

**3/11 STARTING THE METHOD**

Step 1/7 Sieve 1/7

Method name

Coffee Arabica

Sample name

Arabica

Sample ID

06

Comment

Username

GBE

Department

Ok

**4/11 GUIDED SIEVING**  
**STEP 1: SAMPLE IDENTIFICATION**

×

Step 2 / 7 Sieve 1 / 7

Remove sieve  
and lid

Ok

**5/11 GUIDED SIEVING**  
**STEP 2: PREPARATION OF BALANCE**

Step 3/7 Sieve 1/7

Insert the  
sieve - 125 µm

*Please provide a plausible sieve weight  
>=100 g*

Ok

**6/11 GUIDED SIEVING  
STEP 3: WEIGHING THE REQUESTED SIEVE**

Step 4 / 7 Sieve 1 / 7

Fill in the  
sample

Please provide a plausible sample weight

$\geq 0.1$  g

Weighing assistant  
4.63 g - 13.62 g

Ok


**7/11 GUIDED SIEVING**  
**STEP 4: WEIGHING THE SAMPLE. THE WEIGHING ASSISTANT RECOMMENDS A SUITABLE SAMPLE AMOUNT**

Step 5 / 7 Sieve 1 / 7

Put on the lid

*Please provide a plausible lid weight*

*>=200 g*

Ok

**8/11 GUIDED SIEVING  
STEP 5: WEIGHING THE LID**

Step 6 / 7 Sieve 1 / 7

Ready for  
sieving

Duration (mm:ss)

03:00

Pressure (Pa)

2,500 ^  
v

Speed of nozzle drive (rpm)

20.00 ^  
v

Start sieving


**9/11 GUIDED SIEVING**  
**STEP 6: OPTIONAL ADJUSTMENT OF THE SIEVING PARAMETERS AND START**

Step 7 / 7 Sieve 1 / 7

Remove sieve  
and lid

*Please provide a plausible lift-off weight*

<=-531.65 g

Next

**10/11 GUIDED SIEVING**  
**STEP 7: BACKWEIGHING THE RESIDUAL MASS**

Preliminary results

125 µm

Sieving completed successfully

**Results**

Oversize 7.73 g

**Sample information**

Sample amount 10.32 g

Undersize 2.59 g

**Device parameters**

Duration 01:30 mm:ss

Speed of nozzle drive 20 rpm

Pressure 3,000 Pa

Average pressure 2,784 Pa

Pressure std. dev. sigma 877 Pa

Pressure variance 768,927 Pa<sup>2</sup>

Next

11/11 DISPLAY OF RESULTS

AIR JET SIEVING MACHINES AS 200 JET PRO / AS 200 JET PHARMA

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## AIR JET SIEVING MACHINES AS 200 JET PRO / AS 200 JET PHARMA

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
<b>Graphical User Interface, guided sieving</b>	準拠
<b>Mode of operation</b>	Fast sieving, manual sieving and methods
<b>Languages</b>	German and English
<b>Parameter units changeable</b>	Size: mm/µm/mesh/inch Mass: g/kg Time: s/min Pressure: mbar/Pa/psi Volume: ml/cl/l/cm <sup>3</sup>
<b>SOPs</b>	Yes, integrated in methods, more than 40
<b>Individual sample labelling</b>	Yes, for each sample, reference to method during data choice and in overview
<b>Parameter setting per test sieve</b>	準拠
<b>LIMS ready</b>	準拠
<b>possible export format</b>	.xls, .xlsx, .xps, .csv, .pdf, .png, .jpg, .txt, .rtf, html, .mht
<b>Teamviewer online support</b>	準拠
<b>Screenshot function</b>	準拠
<b>Back-up function</b>	準拠
<b>Automatic storage of results</b>	Yes, in machine, auto report function, PDF document or to printer

[www.retsch.jp/AS200jet-pro-pharma](http://www.retsch.jp/AS200jet-pro-pharma)

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### AIR JET SIEVING MACHINES AS 200 JET PRO / JET PHARMA

**Air Jet Sieving Machines for test sieves 203 mm (8") Ø, with automatic vacuum regulation and internal balance, incl. lids for test sieves 1" and 2" height and rubber mallet - please order test sieve and MANDATORY vacuum cleaner (Item No. 22.748.0020 or Item No. 22.748.0021) separately!**

30.034.0001		AS 200 100–240 V, 50/60 Hz jet pro	incl. test report acc. to EN 10204 2.2 evaluation software installed (no EasySieve required)
30.034.0002		AS 200 100–240 V, 50/60 Hz jet pharma (no EasySieve required)	incl. test report acc. to EN 10204 2.2 evaluation software installed incl. AuditTrail & User Management


### VACCUUM CLEANER AS 200 JET PRO / JET PHARMA

22.748.0020	Industrial vacuum cleaner ATTIX 33 230V 50/60Hz (other electrical versions available) MANDATORY vacuum cleaner for 30.034.0001 & 30.034.0002 NOT SUITABLE for 32.935.0008
32.524.0017	HEPA filter for ATTIX 33
32.524.0018	Cover for HEPA filter ATTIX 33
32.524.0019	Replacement filter bags for ATTIX 33, 5 pcs.

### CYCLONE AS 200 JET PRO / JET PHARMA

22.935.0042	High-performance cyclone including holder, collecting receptacle (250 ml and 500 ml sample bottle) and 1.5 m FDA certified hose with clamps (industrial vacuum cleaner 22.748.0020 or 22.748.0021 recommended)
05.707.0208	FDA certified hose, 1.5 m
02.183.0151	Connector for cyclone

### ACCESSORIES AS 200 JET / JET PRO / JET PHARMA

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## ADDITIONAL ITEMS AIR JET SIEVING MACHINES

99.200.0047

IQ/OQ Documentation for AS 200 jet pro

99.200.0046

IQ/OQ Documentation for AS 200 jet pharma

## RECALIBRATION AIR JET SIEVING MACHINES

**The Retsch Air Jet Sieving Machines can be calibrated at Retsch/Germany or at selected subsidiaries/agents**

SER-97.001.0021

Calibration of sieving machines AS 200 jet pro & AS 200 jet pharma  
at Retsch/Germany, incl. test report acc. to EN 10204 2.2