

Obtaining grinding results in the desired size range with reduced fines

By choosing the lowest speed of the cutting mill SM 300 you can reduce the fines in your sample

Cutting mills are the first choice when it comes to size reduction of heterogeneous, soft, medium-hard, fibrous, elastic or tough sample materials such as wood, twigs, electronic scrap or refuse-derived fuels. They are mostly used for preliminary grinding followed by fine grinding, carried out for example in an ultra centrifugal mill like the ZM 200. For some applications, it is desirable to have as much as possible of the sample within a specific size range with the fines being reduced to a minimum. Cutting mills can do the job for bulky samples with the above-mentioned properties. The SM 300 model is particularly suitable due to the variety of accessories (rotors, bottom sieves, etc.) and flexible speed. The latter was extended from 700 min⁻¹ - 3,000 min⁻¹ to 100 min⁻¹ - 3,000 min⁻¹. By selecting the lowest speed when grinding brittle samples which break easily, it is possible to reduce the fines. The low speed also has a positive effect on temperature-sensitive samples.

Correlation between speed and fines

Application example I: 150 g brittle resin pieces with a maximum size of 10 mm were milled in the SM 300 for 2 minutes, using the parallel section rotor and 8 mm bottom sieve. The desired particle size range lay between 2 and 6.3 mm. Starting with 700 min⁻¹, the speed was reduced to 100 min⁻¹ in steps of 200 min⁻¹. Fig. 1 shows that the fine fraction is reduced with decreasing speed. With 700 min⁻¹, only 30% of the sample had a particle size between 2 and 6.3 mm. With 100 min⁻¹, however, 95% of the sample possessed the desired particle size and only 0.5% were too fine (<0.5 mm).




Fig. 1: Correlation between speed of the cutting mill and resulting fine fraction of the sample

Application example II: 130 ml of a brittle NiMo-Al₂O₃ catalyst with initial particle size of 10 mm were ground within 3.5 min by using the 6-disc rotor and various bottom sieves. By selecting appropriate aperture sizes of the sieves and a suitable speed, major part of the sample was reduced to a particle size within the desired range from 0.6 to 1.4 mm (Fig. 2). For the catalyst sample, 80% were obtained in the desired fraction when a 1.5 mm bottom sieve and a speed of 100 min⁻¹ were used (Fig. 3). An additional sieve with smaller aperture sizes was not employed because the sample would then have to be fed very slowly to the hopper to avoid blockages.




Fig. 2: Effect of aperture size of bottom sieve and speed on the grinding result




Fig. 3: Catalyst before and after grinding in the SM 300; the smallest fine fraction was obtained with a 1.5 mm bottom sieve at 100 min⁻¹

Cutting Mill SM 300

The Cutting Mill SM 300 used for these trials is equipped with a 3kW drive with a high torque and RES technology (additional flywheel mass) which allow for particularly effective preliminary size reduction of heterogeneous material mixes. Analytical fineness is usually achieved in one working run. The mill is employed for cutting a large variety of materials. The sample is only slightly warmed during the process, making the SM 300 suitable for grinding heat-sensitive materials. The wide, freely selectable speed range from 100 min⁻¹ to 3000 min⁻¹ helps to substantially reduce the fine fraction. If very light sample material or small amounts are to be processed, the SM 300 can optionally be equipped with the cyclone-suction-combination for optimum results which accepts collecting vessels from 500 ml to 30 l. The wide selection of bottom sieves, hoppers and collecting vessels allows for adaptation to individual application tasks. RETSCH's product portfolio comprises a range of cutting mills, from the economic basic model to the large SM 400 which is suitable for grinding sample pieces with a maximum size of 170 x 220 mm.



Fig 4: Cutting Mill SM 300 (left) and Cutting Mill SM 400 XL (right)