

AS 200 JET PRO & AS 200 JET PHARMA

KEY ADVANTAGES



| Innovative air jet sieving technology

Enables efficient and reproducible **particle size analysis** of fine powders, especially **below 100 µm**.

| All-in-One Device → Sieving → Weighing → Evaluating

Saves lab space and **eliminates errors** from transferring material between devices. **Streamlines workflow** for efficiency.

| Automatic underpressure regulation

Ensures consistent and reliable sieving results by **maintaining optimal air flow** throughout the process.

| User-friendly integrated software and touch display

Easy control of all parameters, fast evaluation, and reporting. Improves user experience with **guided sieving** and process reliability.

| Data export and connectivity

USB and **Ethernet interfaces** for seamless data transfer and integration into **Laboratory Information Management Systems**

| Various checks and assistants

Reduce operating error, ensure reproducible results, and make operation easy even for new users.

| Pharma-conform design

Meets the strict requirements of the pharmaceutical industry, including **full GMP compliance** and easy cleaning thanks to 316L steel sieving chamber. Includes **audit trail, user management** and **e-signatures**.

| IQ/OQ documentation available

Supports qualification and validation processes, **simplifying regulatory compliance**.

SMART.

PRECISE.

EFFICIENT.

